

# Experience treatment with **NO MORE NEEDLES**

If you are looking to be free from needles, **ZOMA-Jet 10\*** may be right for you

## Reusable needle-free head

Each needle-free head can be used for up to 7 days

## Easy viewing dosage window

Displays dosage for accurate ZOMACTON® (somatropin) for Injection 10mg dose delivery

## One-push dose-delivery button

Go from needle-sticks to quick clicks

Absolutely  
**NO  
NEEDLES**



Actual size

## Needle-free ZOMA-Jet 10 also:

- **Eliminates concerns** about needle-sticks and fear of needles
- **Involves 3 routine steps** for completely needle-free dose delivery after you become familiar with the instructions for use
- **Requires 1 quick click** for dose delivery
- **Offers initial in-home training** by nurse educators to help you get started with therapy

NEEDLE-FREE

**ZOMA-Jet 10**

**NO STICK,**

**QUICK CLICK**



To learn how to get started with **Needle-free ZOMA-Jet 10**, speak with a ZOGO Care Coordinator by calling **844-ZomaJet (966-2538)** (Mon-Fri, 9AM to 7PM EST)






\*For Zomacton 10mg only.

# Choose the device that meets your specific needs

There are a number of devices available. Choosing the right one is important.

## FEATURES OF AVAILABLE DEVICE OPTIONS

● =Yes    EMPTY BOX=No

Product and Device	Number of Indications	Prefilled	Product Storage*	Product Refrigeration Required*	Priming Required	Dial Backwards	Needle-free Dose Delivery	Hold Time
<b>ZOMACTON®</b> <b>ZOMA-Jet® 10<sup>1</sup></b> 	6		Up to 28 days	●	●	●	●	2 seconds
<b>Humatrope®</b> <b>HumatroPen®2</b> 	6		Up to 28 days	●	●	●		5 seconds
<b>Norditropin®</b> <b>FlexPro®3</b> 	7	●	Up to 21 days†		●	●		6 seconds
<b>Omnitrope®</b> <b>pen<sup>4</sup></b> 	6		Up to 28 days	●	●			5 seconds
<b>Nutropin AQ®</b> <b>NuSpin®5</b> 	5	●	Up to 28 days	●	●	●		5 seconds
<b>Genotropin®</b> <b>Pen®6</b> 	6		Up to 28 days	●	●	●		5 seconds
<b>Genotropin</b> <b>MiniQuick®6</b> 	6	●	Up to 24 hours	●				A few seconds
<b>Saizen</b> <b>one.click®7,8</b> 	2		Up to 21 days	●		●		10 seconds
<b>Saizen</b> <b>easypod®7,9</b> 	2		Up to 21 days	●				5 to 10 seconds

\*After reconstitution/initial use.

†Up to 28 days if refrigerated.

The information set forth herein is presented to provide reference to a non-exhaustive list of various characteristics of the devices listed above. This information does not suggest or imply in any manner that any of the characteristics of each of the devices, alone, or in combination, are superior over another, nor is it intended to imply that the products or devices mentioned have comparable safety or efficacy, indications, and doses. Indications are not the same across products.

**References:** 1. ZOMACTON [prescribing information], Parsippany, NJ: Ferring Pharmaceuticals. 2. Humatrope [prescribing information], Indianapolis, IN: Eli Lilly and Company. 3. Norditropin [prescribing information], Plainsboro, NJ: Novo Nordisk A/S. 4. Omnitrope [prescribing information], Princeton, NJ: Sandoz, A Novartis Division. 5. Nutropin AQ [prescribing information], San Francisco, CA: Genentech, Inc. 6. Genotropin [prescribing information], New York, NY: Pharmacia & Upjohn Co, Division of Pfizer Inc. 7. Saizen [prescribing information], Rockland, MA: EMD Serono, Inc. 8. Saizen one.click [instructions for use], Rockland, MA: EMD Serono, Inc. 9. Saizen easypod [instructions for use], Rockland, MA: EMD Serono, Inc.

Genotropin, Genotropin Pen, and Genotropin MiniQuick are registered trademarks of Pfizer; Humatrope and HumatroPen are registered trademarks of Eli Lilly and Company; Norditropin and FlexPro are registered trademarks of Novo Nordisk Health Care AG; NuSpin and Nutropin AQ are registered trademarks of Genentech, Inc; Omnitrope is a registered trademark of Novartis AG; Saizen, Saizen one.click, and Saizen easypod are registered trademarks of EMD Serono, Inc.